


Basics of Color Doppler Ultrasound and Applications

Introduction

Color Doppler is the most demanding diagnostic technology for doctors in the current scenario. It gives the information of flowing blood in terms of its Velocity, Direction, Acceleration Time duration etc. Generally it indicates the blood flow direction with Red color if the blood is flowing towards the transducer, and with Blue color when the blood is moving away from the transducer. Doctors are able to diagnose many diseases with the help of Color Doppler technology.

Basic Principle

Doppler ultrasound measures the movement of the scatterers through the beam as phase change in the received signal. The resulting Doppler frequency can be used to measure velocity if the beam/flow angle is known.



The size of the Doppler signal is dependent on:

- 1. Blood velocity
- 2. Ultrasound frequency

3. Angle of incidence




TEE (Trans Esophagus Echocariography) probe



Applications: cardiac

Liner probe

Applications: Breast, small parts, nerve, vascular, venous access, musculoskeletal,



Types of Transducers

Convex probe



Applications: Obstetrics, gynecology

Endorectal probe

Applications: Obstetrics, gynecology

(A) Higher-frequency Doppler signal is obtained if the beam is aligned more to the direction of flow. In the diagram, beam (A) is more aligned than (B) and produces higher-frequency Doppler signals. The beam/flow angle at (C) is almost 90° and there is a very poor Doppler signal. The flow at (D) is away from the beam and there is a negative signal.

Types of Doppler Color Doppler (CD)

conditions.

Color Power Doppler (CPD) CPD is a type of Color Doppler to visualize the presence of detectable blood flow. The flow information is based on the amplitude or strength of echoes received from moving cells and not frequency shifts.

Directional Color Power Doppler (DCPD) Directional Color Power Doppler combines power (amplitude) of Doppler signal with directional (phase) information to encode direction and variations in blood flow.

Spectral Doppler :

Spectral Doppler refers to the combination of either continuous wave Doppler or pulsed Doppler with a spectral display. Spectral Doppler provides a quantitative analysis of the velocity and direction of blood flow.

Continuous Wave Doppler

Continuous Wave (CW) Doppler is an ultrasound imaging mode, which records blood flow velocities along the length of the beam. Continuous wave Doppler uses different crystals to send and receive the signal.

Pulse Wave Doppler (PWD): Evaluates blood flow velocities in a range specific area along the length of the sound beam.

Transcranial Doppler (TCD): Transcranial Color Doppler sonography allows evaluating the presence and flow direction of vessels as well as their relationships to surrounding structures.

Imaging Modes:

A (Amplitude) Mode For Eyes B/BB (Bright / Double Bright) Mode For Abdominal & Pelvic organs BC (Compound Bright) Mode For Abdominal & Pelvic organs M (Motion) Mode

Advance Imaging	
DPDBF	(Do
THI	(Tis
TDI	(Tis
TSS	(Tis
SRA	(Spe
1	


Different Medical Applications –

Ultrasound can be applicable for different medical analysis like, Abdominal, Obstetrics, Gynecology, Urology, Cardiology, Vascular, Musculoskeletal and Small body parts











Ventricular septum TDI





IT-2013, Ramchandrapura Indl. Area, Sitapura Extension, Jaipur-302022, India. Ph: +91-9829132777; +91-9413330765; Email: info@tesca.in; Web: www.tesca.in


Color Doppler is an ultrasound imaging mode, which visualizes the presence, direction and velocity of flowing blood in a wide range of flow

For Cardiac Analysis

g Processing Technology:

uble Phase Digital Beam Forming) sue Harmonics Imaging) sue Doppler Imaging) sue Speed of Sound) eckle Reduction Algorithm)

LV short axis



71800 **Bio Medical Instrumentation**

