

### Features

- Compact, mobile unit, designed to provide required services for testing various Heat Exchangers using single source.
- Entire process is computer controlled, including change of direction of flow for parallel & counter-flow analysis.
- Comprehensive Instrumentation Panel with all necessary safety instruments.

Heat exchangers transfer thermal energy from the flow of one medium to another. The two flows do not come into direct contact with one another. Efficient heat transfer is a prerequisite for economical processes. Therefore, different heat exchanger types are used in practice depending on the requirements.

Tesca Computer Controlled Heat Exchanger Service 32346 can be used to investigate and compare different heat exchanger designs. The complete experimental setup consists of two main elements: 32346 as supply and control unit and choice of heat exchanger: Tubular heat exchanger (32347), plate heat exchanger (32348), shell and tube heat exchanger (32350) and stirred tank with jacketed vessel and coil (32351), and Linear Heat Transfer Unit 32352. Water is used as the medium.

The heat exchanger to be investigated is connected to the supply unit. The hot water flows through the heat exchanger. Part of the thermal energy of the hot water is transferred to the cold water.

Reversing the water connections changes the direction of flow and thus allows parallel flow or counter-flow operation.

The main function of the 32346 is to provide the required cold and hot water circuits. To do this, the supply unit is equipped with a heated tank and pump for the hot water circuit, connections for the cold water circuit and a switch cabinet with displays and controls. A temperature controller controls the hot water temperature. The flow rate in the hot water and cold water circuit is adjusted using valves. The cold water circuit can be fed from the laboratory mains or from Water Chiller.

Optionally, Tesca offers software for data acquisition and an educational software. With explanatory texts and illustrations, the educational software significantly aids the understanding of the theoretical principles. With the aid of an authoring system, the teacher can create further exercises.

Sensors record the temperatures and flow rates. The measured values are read from digital displays and can be transmitted simultaneously via USB directly to a PC where they can be analyzed using the software.

### Specifications

- Supply unit for heat exchangers
- Hot water circuit with tank, heater, temperature

Note: Specifications are subject to change.

### Tesca Technologies Pvt. Ltd.

IT-2013, Ramchandrapura Industrial Area, Sitapura Extension,  
Near Bombay Hospital, Vidhani Circle, Jaipur-302022, Rajasthan, India,  
Tel: +91-141-2771791 / 2771792; Email: info@tesca.in, tesca.technologies@gmail.com  
Website: www.tescaglobal.com



controller, pump and protection against lack of water

- Cold water circuit from laboratory mains or water chiller
- Temperature controller controls the temperature of hot water
- Flow adjustable using valves
- Digital displays for 6 temperature and 2 flow rate sensors
- Water connections with quick-release couplings
- Stirring machine connection with speed adjustment (WL 110.04)
- Optional Tesca software: educational software and data acquisition
- Software for data acquisition via USB under Windows 7, 8.1, 10

### Technical Specifications

- Pump
  - power consumption: 120W
  - max. flow rate: 600L/h
  - max. head: 30m
- Heater
  - Power output: 3kW
  - Thermostat: 0...70°C
  - Hot water tank: approx. 10L
- Measuring ranges
  - Temperature: 6x 0...100°C
  - Flow rate: 2x 20...250L/h

### Interface In-built Module (Optional):

This control interface is common for the 'Tesca' Heat Exchangers and can work with one or several exchangers.

The Control Interface is part of the SCADA system. Control interface with process diagram on the front panel.

The unit control elements are can be computer controlled.

- Simultaneous visualization in the computer of all parameters involved in the process.
- Calibration of all sensors involved in the process.
- Real time curves representation about system responses.
- All the actuators' values can be changed at any time from the keyboard allowing the analysis about curves and responses of the whole process.
- Shield and filtered signals to avoid external interferences.
- Real time PID control with flexibility of modifications from the computer keyboard of the PID parameters, at any moment during the process.
- Real time PID control for parameters involved in the process simultaneously.
- Proportional control, integral control and derivative control, based on the real PID mathematical formula, by changing the values, at any time, of the three control constants (proportional, integral and derivative constants).
- Open control allowing modifications, at any moment and in real time, of parameters involved in the process simultaneously.

Three safety levels, one mechanical in the unit, another electronic in the control interface and the third one in the control software.

### OPTIONAL HEAT EXCHANGER MODULES

#### Tubular Heat Exchanger Module

Order Code - 32347



#### Specifications

- Tubular heat exchanger for connection to 32346
- Hot and cold water supply from 32346
- Parallel flow and counter-flow operation possible
- Recording of temperature using 32346 and two additional temperature sensors for measuring the central temperature

#### Technical Specifications

- Heat transfer surfaces
  - Mean transfer surface: 250cm<sup>2</sup>
- Inner tube, stainless steel

Note: Specifications are subject to change.

### Tesca Technologies Pvt. Ltd.

IT-2013, Ramchandrapura Industrial Area, Sitapura Extension,  
Near Bombay Hospital, Vidhani Circle, Jaipur-302022, Rajasthan, India,  
Tel: +91-141-2771791 / 2771792; Email: info@tesca.in, tesca.technologies@gmail.com  
Website: www.tescaglobal.com

- Outer diameter: 12mm
  - Wall thickness: 1mm
  - Outer tube, transparent (PMMA)
    - Outer diameter: 20mm
    - Wall thickness: 2mm
  - Measuring ranges
    - Temperature: 2x 0...100°C
- Easy connection to the Base Service Unit 32346.

#### Optional


Analogue and digital PID control. PID menu and set point selection required in the whole work range.

- Management, processing, comparison and storage of data.
- Sampling velocity up to 250 KS/s (kilo samples per second).
- Calibration system for the sensors involved in the process.
- Allows the registration of the alarms state and the graphic representation in real time.
- Comparative analysis of the obtained data, after the process and modification of the conditions during the process.
- Software part of the SCADA system from 'Interface In-built Module'

#### Plate Heat Exchanger Module

Order Code - 32348

Tesca Plate Heat Exchanger Module 32348 allows the study of heat transfer between hot and cold water through alternate channels formed between parallel plates. The exchanger allows measuring cold and hot temperatures at the inlet and outlet of the exchanger.



#### Specifications

- Plate heat exchanger for connection to 32346
- Hot and cold water supply from 32346
- Parallel flow and counter-flow operation possible
- Six soldered plates
- recording of temperature using 32346

#### Technical Specifications

- Four ports or connections of hot and cold water input and output.
- Maximum flow: 12m<sup>3</sup>/h.
- Maximum work pressure: 10 bar.
- Maximum work temperature: 100 °C.
- Minimum work temperature: 0 °C.
- Maximum number of plates: 20.
- Internal circuit capacity: 0.176 l.

- External circuit capacity: 0.22 l.
- Area: 0.32 m<sup>2</sup>.
- Four temperature sensors ("J" type):
  - Two temperature sensors for measuring cold water temperature (inlet and outlet).
  - Two temperature sensors for measuring hot water temperature (inlet and outlet).
- Easy connection to the Base Service Unit 32346.

#### Optional

Analog and digital PID control. PID menu and set point selection required in the whole work range.

- Management, processing, comparison and storage of data.
- Sampling velocity up to 250 KS/s (kilo samples per second).
- Calibration system for the sensors involved in the process.
- Allows the registration of the alarms state and the graphic representation in real time.
- Comparative analysis of the obtained data, after the process and modification of the conditions during the process.
- Software part of the SCADA system from 'Interface In-built Module'

#### Extended Plate Heat Exchanger Module Order Code - 32349

Tesca Extended Plate Heat Exchanger Module 32349 allows the study of heat transfer between hot and cold water through alternate channels formed between parallel plates.

The exchanger allows measuring cold and hot temperatures at different points of the exchanger.

#### Technical Specifications

- Four ports or connections of hot and cold water input and output.
- Maximum flow: 12 m<sup>3</sup>/h.
- Maximum work pressure: 10 bar.
- Maximum work temperature: 100 °C.
- Minimum work temperature: 0 °C.
- Maximum number of plates: 20.
- Internal circuit capacity: 0.176 l.
- External circuit capacity: 0.22 l.
- Area: 0.32 m<sup>2</sup>.
- Ten temperature sensors ("J" type):
  - Five temperature sensors for measuring cold water temperature (inlet, outlet and interim positions).
  - Five temperature sensors for measuring hot water temperature (inlet, outlet and interim positions).
- Easy connection to the Base Service Unit 32346.

#### Optional

- Analog and digital PID control. PID menu and set point selection required in the whole work range.

Note: Specifications are subject to change.

### Tesca Technologies Pvt. Ltd.

IT-2013, Ramchandrapura Industrial Area, Sitapura Extension,  
Near Bombay Hospital, Vidhani Circle, Jaipur-302022, Rajasthan, India,  
Tel: +91-141-2771791 / 2771792; Email: info@tesca.in, tesca.technologies@gmail.com  
Website: www.tescaglobal.com

- Management, processing, comparison and storage of data.
- Sampling velocity up to 250 KS/s (kilo samples per second).
- Calibration system for the sensors involved in the process.
- Allows the registration of the alarms state and the graphic representation in real time.
- Comparative analysis of the obtained data, after the process and modification of the conditions during the process.
- Software part of the SCADA system from 'Interface In-built Module'

#### Shell and Tube Heat Exchanger Module Order Code - 32350

Tesca Shell and Tube Heat Exchanger Module 32350 consists on a group of tubes inside the heat exchanger. The hot water flows through the internal tubes and the cooling water circulates through the space between the internal tubes and the shell.



There are traverse baffles placed in the shell to guide the cold water to maximize the heat transfer.

#### Specifications

- Shell and tube heat exchanger (cross-flow) for connection to 32346
- Hot and cold water supply from 32346
- Cross parallel flow and cross counter-flow operation possible
- Transparent shell, visible tube bundle
- Tube bundle consisting of 7 tubes and 4 baffle plates
- Recording of temperature using 32346

#### Technical Specifications

Four segmented baffles located transversally in the shell.

- Exchange length of the shell and each tube: L = 0.5 m.
- Internal tube (21 tubes):
  - Internal diameter: D<sub>int</sub> = 8 • 10<sup>-3</sup> m.
  - External diameter: D<sub>ext</sub> = 10 • 10<sup>-3</sup> m.
  - Thickness = 10<sup>-3</sup> m.
  - Internal heat transfer area: A<sub>h</sub> = 0.0126 m<sup>2</sup>.
  - External heat transfer area: A<sub>c</sub> = 0.0157 m<sup>2</sup>.
- Shell:
  - Internal diameter: D<sub>int,c</sub> = 0.148 m.
  - External diameter: D<sub>ext,c</sub> = 0.160 m.
  - Thickness = 6 • 10<sup>-3</sup> m.
- Seven temperature sensors ("J" type), for measuring cold and hot water temperatures at

- different points of the exchanger.
- Easy connection to the Base Service Unit 32346.

#### Optional:

- Analog and digital PID control. PID menu and set point selection required in the whole work range.
- Management, processing, comparison and storage of data.
  - Sampling velocity up to 250 KS/s (kilo samples per second).
  - Calibration system for the sensors involved in the process.
  - Allows the registration of the alarms state and the graphic representation in real time.
  - Comparative analysis of the obtained data, after the process and modification of the conditions during the process.
  - Software part of the SCADA system from 'Interface In-built Module'

#### Jacketed Vessel Heat Exchanger Module Order Code - 32351



Tesca Jacketed Vessel Heat Exchanger Module 32351 allows the study of heat transfer between hot water flowing through a jacket and the cold water contained in a vessel.

It can work in continuous supply or in a batch process (heating of a constant mass of water contained in a vessel).

The exchanger allows measuring temperatures at the inlet and outlet of the exchanger in cold as well as in hot water.

#### Specifications

- Stirred tank for connection to 32346
- Hot and cold water supply from 32346
- Heating using jacket or coiled tube
- Stirring machine can be used in all modes
- Speed of stirring machine adjustable using 32346
- Visible working area due to transparent cover
- Recording of temperature using 32346 and additional temperature sensor for measuring temperature in tank

#### Technical Specifications

- Main metallic elements made of stainless steel.

Note: Specifications are subject to change.

#### Tesca Technologies Pvt. Ltd.

IT-2013, Ramchandrapura Industrial Area, Sitapura Extension,  
Near Bombay Hospital, Vidhani Circle, Jaipur-302022, Rajasthan, India,  
Tel: +91-141-2771791 / 2771792; Email: info@tesca.in, tesca.technologies@gmail.com  
Website: www.tescaglobal.com

- Diagram in the front panel with distribution of the elements similar to the real one.
- Constituted of a vessel.
- Vessel total volume: 14 l.
- Internal vessel volume: 7 l approx.
- Jacket volume: 7 l approx.
- Overflow or a pipe allows the exit of the water in the vessel through its upper part to maintain a constant flow during the process with a continuous supply.
- Jacket surrounds the vessel through where hot water flows. Electric stirrer.
- Five temperature sensors ("J" type):
  - Three temperature sensors for measuring cold water temperature.
  - Two temperature sensors for measuring hot water temperature.
- Easy connection to the Base Service Unit 32346.

#### Optional

Analog and digital PID control. PID menu and set point selection required in the whole work range.

- Management, processing, comparison and storage of data.
- Sampling velocity up to 250 KS/s (kilo samples per second).
- Calibration system for the sensors involved in the process.
- Allows the registration of the alarms state and the graphic representation in real time.
- Comparative analysis of the obtained data, after the process and modification of the conditions during the process.
- Software part of the SCADA system from 'Interface In-built Module'

#### Experiments

- Demonstration of indirect heating/ cooling by transfer of heat from one fluid stream to another when separated by a solid wall.
- Energy balance determination (heat balance) and calculation of efficiencies by measuring the flow rates and temperature changes in the hot and cold fluid streams.
- Introduction to different styles of heat exchanger and comparison of the differences in operation and performance.
- Using the Logarithmic Mean Temperature Difference (LMTD) in heat transfer calculations.
- Definition and measurement of Overall Heat Transfer Coefficient (U).
- Demonstration of the differences between counter-current and co-current operation.
- Demonstration of the transition from linear to turbulent flow.
- Effect of hot and cold fluid flow rate on the heat transfer coefficient.
- Effect of driving force (temperature differential)

- on the heat transfer coefficient.
- Investigation of heat loss and reduction in heat transfer coefficient due to fouling of the heat transfer surfaces (suitable student project using user induced. fouling).

### Linear Conduction Heat Transfer Order Code - 32352



Tesca Linear Conduction Heat Transfer 32352 has been designed for students to study the phenomena of Linear Heat Conduction. The setup consists of a Heating Section, Cooling Section & Specimens of different metals. The specimens can be fitted in between Heating & Cooling section. Test specimens of same metal or different metals can be clamped in between the heating & cooling sections. Temperature sensors record the surface temperature along length of the specimen. Instruments are provided to measure the Temperatures, Power Input to heater & Cooling Water Flow rate. Detailed Operation & Maintenance Manual is provided along with the trainer.

#### Specifications

- Power Supply with power regulator
- Heating Section: MOC Brass, 25 mm Dia, fitted with Cartridge type Heater Capacity – 100 W.
- Cooling section: MOC Brass, 25 mm Dia, with cooling water jacket.
- Test Specimen:
  - Brass section of 25 mm diameter and fitted with two thermocouples at the same intervals.
  - Stainless steel section of 25 mm diameter.
  - Aluminium section of 25 mm diameter
  - Brass section reduced in diameter to 13mm
- Steel Frame
- Thermocouple sensors
- Temperature Indicator
- Variable Area water flow meter
- Digital Voltmeter
- Digital Ammeter
- Easy connection to the Base Service Unit 32346.

#### Optional

- Analog and digital PID control. PID menu and set

Note: Specifications are subject to change.

### Tesca Technologies Pvt. Ltd.

IT-2013, Ramchandrapura Industrial Area, Sitapura Extension,  
Near Bombay Hospital, Vidhani Circle, Jaipur-302022, Rajasthan, India,  
Tel: +91-141-2771791 / 2771792; Email: info@tesca.in, tesca.technologies@gmail.com  
Website: www.tescaglobal.com

- point selection required in the whole work range.
- Management, processing, comparison and storage of data.
- Sampling velocity up to 250 KS/s (kilo samples per second).
- Calibration system for the sensors involved in the process.
- Allows the registration of the alarms state and the graphic representation in real time.
- Comparative analysis of the obtained data, after the process and modification of the conditions during the process.
- Software part of the SCADA system from 'Interface In-built Module'

#### Experiment Capabilities

- Understanding the Fourier rate equation in determining the rate of heat flow through solid materials
- Measuring the temperature distribution for steady state conduction of energy through a uniform plane solid and a composite plane solid
- Determine the constant of proportionality (thermal conductivity) of different materials (conductors and insulators)
- Measuring the temperature drop at the contact face between adjacent layers in a composite plane solid
- Measuring the temperature distribution for steady state conduction of energy through a plane solid of reduced cross sectional area
- Understanding the application of poor conductors (insulators)
- Observing unsteady state conduction (qualitative only)

### Concentric Tube Heat Exchanger Module Order Code - 32353



Tesca Concentric Tube Heat Exchanger Module allows the study of heat transfer between hot water flowing through an internal tube and cold water flowing in the ring area lying between the internal and the external tubes.

This exchanger allows measuring hot and cold water temperatures at different points of the exchanger.

#### Technical Specifications

- Exchange length:  $L = 2 \times 0.5 = 1$  m.
- Internal tube:
- Internal diameter:  $D_{int} = 16 \cdot 10^{-3}$  m.

- External diameter:  $D_{ext} = 18 \cdot 10^{-3}$  m.
- Thickness =  $10^{-3}$  m.
- Heat transfer internal area:  $A_h = 0.0503$  m<sup>2</sup>.
- Heat transfer external area:  $A_c = 0.0565$  m<sup>2</sup>.
- External tube:
- Internal diameter:  $D_{int} = 26 \cdot 10^{-3}$  m.
- External diameter:  $D_{ext} = 28 \cdot 10^{-3}$  m.
- Thickness =  $10^{-3}$  m.
- Six temperature sensors ("J" type)
- Extended Concentric Tube Heat Exchanger:
- Exchange length:  $L = 4 \times 1 = 4$  m.
- Internal tube:
- Internal diameter:  $D_{int} = 16 \cdot 10^{-3}$  m.
- External diameter:  $D_{ext} = 18 \cdot 10^{-3}$  m.
- Thickness =  $10^{-3}$  m.
- Heat transfer internal area:  $A_h = 0.0503$  m<sup>2</sup>.
- Heat transfer external area:  $A_c = 0.0565$  m<sup>2</sup>.
- External tube:
- Internal diameter:  $D_{int} = 26 \cdot 10^{-3}$  m.
- External diameter:  $D_{ext} = 28 \cdot 10^{-3}$  m.
- Thickness =  $10^{-3}$  m.
- Ten temperature sensors ("J" type):
  - Three temperature sensors for measuring cold water temperature:
    - Cold water inlet.
    - Cold water mid-position.
    - Cold water outlet.
  - Three temperature sensors for measuring hot water temperature:
    - Hot water inlet.
    - Hot water mid-position.
    - Hot water outlet.
- Easy connection to the Base Service Unit 32346.

### Optional

Analog and digital PID control. PID menu and set point selection required in the whole work range.

- Management, processing, comparison and storage of data.
- Sampling velocity up to 250 KS/s (kilo samples per second).
- Calibration system for the sensors involved in the process.
- Allows the registration of the alarms state and the graphic representation in real time.
- Comparative analysis of the obtained data, after the process and modification of the conditions during the process.
- Software part of the SCADA system from 'Interface In-built Module'

### Extended Concentric Tube Heat Exchanger Module : Order Code - 32354

Tesca Extended Concentric Tube Heat Exchanger Module 32354 allows the study of heat transfer between hot water flowing through an internal tube and cold water flowing in the ring area lying between the internal and the external tubes.

Note: Specifications are subject to change.

### Tesca Technologies Pvt. Ltd.

IT-2013, Ramchandrapura Industrial Area, Sitapura Extension,  
Near Bombay Hospital, Vidhani Circle, Jaipur-302022, Rajasthan, India,  
Tel: +91-141-2771791 / 2771792; Email: info@tesca.in, tesca.technologies@gmail.com  
Website: www.tescaglobal.com

This exchanger allows measuring hot and cold water temperatures at different points of the exchanger. Is a more sophisticated unit than 32353, with four longer tube sections, giving four times the overall heat transfer area and three interim temperature measurement points (temperature sensors) in each fluid stream.

This exchanger has sufficient heat transfer area for demonstrating the typical counter-current flow conditions where the outlet of the heated stream is hotter than the outlet of the cooled stream.

The exchanger is formed by two concentric copper tubes with hot water circulating through the interior tube and cold water circulating in the ring space.

### Technical Specifications

- Four equal sections of 1000 mm each one, where heat transfer takes place.
- Exchange length: 1 mtr x 4 nos.
- Internal tube:
  - Internal diameter:  $D_{int} = 16 \cdot 10^{-3}$  m.
  - External diameter:  $D_{ext} = 18 \cdot 10^{-3}$  m.
  - Thickness =  $10^{-3}$  m.
  - Heat transfer internal area:  $A_h = 0.0503$  m<sup>2</sup>.
  - Heat transfer external area:  $A_c = 0.0565$  m<sup>2</sup>.
- External tube:
  - Internal diameter:  $D_{int} = 26 \cdot 10^{-3}$  m.
  - External diameter:  $D_{ext} = 28 \cdot 10^{-3}$  m.
  - Thickness =  $10^{-3}$  m.
- Ten temperature sensors ("J" type):
  - \* Five temperature sensors for measuring cold water temperature:
    - Cold water inlet.
    - Cold water at different interim positions (3).
    - Cold water outlet.
  - \* Five temperature sensors for measuring hot water temperature:
    - Hot water inlet.
    - Hot water at different interim positions (3).
    - Hot water outlet.
- Easy connection to the Base Service Unit 32346.

### Optional

Analog and digital PID control. PID menu and set point selection required in the whole work range.

- Management, processing, comparison and storage of data.
- Sampling velocity up to 250 KS/s (kilo samples per second).
- Calibration system for the sensors involved in the process.
- Allows the registration of the alarms state and the graphic representation in real time.
- Comparative analysis of the obtained data, after the process and modification of the conditions during the process.
- Software part of the SCADA system from 'Interface In-built Module'

### Coil Vessel Heat Exchanger Module

#### Order Code - 32355

Tesca Coil Vessel Heat Exchanger Module 32355 allows the study of heat transfer between hot water flowing through a coil and cold water contained in the vessel.

It can work in continuous supply or in a batch process.

#### Technical Specifications

- PVC -glass vessel, volume: 14 l.
- Overflow or pvc-glass tube lets the output of water from the vessel in the upper part in order to maintain the flow constant for continue supply process.
- Copper coil where the water circulates:
  - Dint= 4.35 mm.
  - Dext= 6.35 mm.
  - Total length of the tube that forms the coil: 1.5 m.
- Electric stirrer.
- Five temperature sensors ("J" type):
  - Three temperature sensors for measuring cold water temperature.
  - Two temperature sensors for measuring hot water temperature.
- Easy connection to the Base Service Unit 32346.

#### Optional

Analog and digital PID control. PID menu and set point selection required in the whole work range.

- Management, processing, comparison and storage of data.
- Sampling velocity up to 250 KS/s (kilo samples per second).
- Calibration system for the sensors involved in the process.
- Allows the registration of the alarms state and the graphic representation in real time.
- Comparative analysis of the obtained data, after the process and modification of the conditions during the process.
- Software part of the SCADA system from 'Interface In-built Module'

### Turbulent Flow Heat Exchanger Module

#### Order Code - 32356

Tesca Turbulent Flow Heat Exchanger Module 32356 lets us study the heat transfer between hot water that circulates through an internal tube and cold water that flows through the annular zone between the internal and the external tube. This exchanger let us measure cold water and hot water temperatures at different points of the exchanger.

#### Technical Specifications

- Formed by two copper concentric tubes with hot water circulating through the internal tube and cold water circulating through the annular space.

Note: Specifications are subject to change.

### Tesca Technologies Pvt. Ltd.

IT-2013, Ramchandrapura Industrial Area, Sitapura Extension,  
Near Bombay Hospital, Vidhani Circle, Jaipur-302022, Rajasthan, India,  
Tel: +91-141-2771791 / 2771792; Email: info@tesca.in, tesca.technologies@gmail.com  
Website: www.tescaglobal.com

- The exchanger has 4 equal sections of 500 mm each one, where the heat transfer takes place.
- Exchange length:  $L = 4 \times 0.5 = 2$  m.
- Internal tube:
  - Internal diameter:  $D_{int} = 8 \cdot 10^{-3}$  m.
  - External diameter:  $D_{ext} = 10 \cdot 10^{-3}$  m.
  - Thickness =  $10^{-3}$  m.
  - Internal heat transfer area:  $A_h = 0.0377$  m<sup>2</sup>.
  - External heat transfer area:  $A_c = 0.0471$  m<sup>2</sup>.
- External tube:
  - Internal diameter:  $D_{int,c} = 13 \cdot 10^{-3}$  m.
  - External diameter:  $D_{ext,c} = 15 \cdot 10^{-3}$  m.
- Thickness =  $10^{-3}$  m.
- Twelve temperature sensors ("J" type):
  - Cold water temperature sensor at the exchanger inlet or outlet.
  - Hot water sensor at the exchanger inlet.
  - Cold water sensor between the first and second stretch of the exchanger.
  - Hot water sensor between the first and second stretch of the exchanger.
  - Cold water sensor between the second and third stretch of the exchanger.
  - Hot water sensor between the second and third stretch of the exchanger.
  - Cold water sensor between the third and fourth stretch of the exchanger.
  - Hot water sensor between the third and fourth stretch of the exchanger.
  - Cold water temperature sensor at the exchanger inlet or outlet.
  - Hot water sensor at the exchanger outlet.
  - Temperature sensor of the exterior surface of the interior tube at the exchanger inlet.
  - Temperature sensor of the exterior surface of the interior tube at the exchanger outlet.
- Easy connection to the Base Service Unit 32346.

#### Optional

Analog and digital PID control. PID menu and set point selection required in the whole work range.

- Management, processing, comparison and storage of data.
- Sampling velocity up to 250 KS/s (kilo samples per second).
- Calibration system for the sensors involved in the process.
- Allows the registration of the alarms state and the graphic representation in real time.
- Comparative analysis of the obtained data, after the process and modification of the conditions during the process.
- Software part of the SCADA system from 'Interface In-built Module'

### **Cross Flow Heat Exchanger Module Order Code - 32357**

Tesca Cross Flow Heat Exchanger Module 32357 is designed to study heat transfer between two fluids in cross flow configuration.

Hot water flow coming from the base unit enters and leaves a radiator perpendicular to an air current, which is generated by a fan.

The heat exchanger allows to measure water and air temperatures at the inlet and outlet of the exchanger.

#### **Technical Specifications**

- PMMA rectangular duct of 800 x 200 x 200 mm.
- Radiator located across the air duct.
- The fins of the radiator are made of aluminum and have a heat transfer area of 35000 mm<sup>2</sup>.
- Axial fan with speed control from computer (PC). It provides a maximum air velocity of 3 m/s.
- Four "J" type temperature sensors to measure input and output water and air temperatures.
- Velocity sensor to measure air velocity, range: 0 – 4 m/s.
- Two ball valves.
- Easy connection to the Base Service Unit 32346.

#### **Optional**

Analog and digital PID control. PID menu and set point selection required in the whole work range.

- Management, processing, comparison and storage of data.
- Sampling velocity up to 250 KS/s (kilo samples per second).
- Calibration system for the sensors involved in the process.
- Allows the registration of the alarms state and the graphic representation in real time.
- Comparative analysis of the obtained data, after the process and modification of the conditions during the process.
- Software part of the SCADA system from 'Interface In-built Module'

#### **Accessories & Manuals**

- All trainers & modules are provided with cables and accessories for normal operation.
- Manuals to include Services, Assembly and Installation, Interface and Control Software, Starting-up, Safety, Maintenance, Calibration & Practices are provided.

#### **Services Required**

- Electric Supply 230 V AC, 16 A, Single Phase, Earthed.
- Cold Water Supply & Drainage.
- Personal Computer with USB port, Windows operating system & all peripherals.

Note: Specifications are subject to change.

### **Tesca Technologies Pvt. Ltd.**

IT-2013, Ramchandrapura Industrial Area, Sitapura Extension,  
Near Bombay Hospital, Vidhani Circle, Jaipur-302022, Rajasthan, India,  
Tel: +91-141-2771791 / 2771792; Email: info@tesca.in, tesca.technologies@gmail.com  
Website: www.tescaglobal.com

