



### Features

- Demonstration of fundamental aspects of mechanical vibration
- Damping and resonance with forced vibrations
- Two different principles of exciting vibration

Tesca Free & Forced Vibration Apparatus is an experimental unit based on the theory that vibrations are rapidly and securely assembled on a double profile frame made of aluminum using slot nuts and clamping levers. The oscillator consists of a beam mounted in ball bearings at one end; a helical spring is hung on the other end of the beam. The attachment of springs, the exciter and a damper to a perforated panel permits a wide range of different set-ups to be reproduced. Either an unbalance exciter or a displacement exciter generate vibrations; the frequency of the vibrations can be adjusted using the electronic control unit. The

Note: Specifications are subject to change.

### **Tesca Technologies Pvt. Ltd.**

IT-2013, Ramchandrapura Industrial Area, Sitapura Extension,  
Near Bombay Hospital, Vidhani Circle, Jaipur-302022, Rajasthan, India,  
Tel: +91-141-2771791 / 2771792; Email: info@tesca.in, tesca.technologies@gmail.com  
Website: www.tescaglobal.com

displacement exciter can be fitted directly to the base of the spring. The vibrations can be damped using the adjustable viscosity damper. To record vibration processes over time, a mechanical drum plotter is included. The experiment also includes an amplitude contact with a TTL output, e.g. for triggering stroboscopes.

The measured values can be evaluated with the software for data acquisition (optional). A practice set for torsional vibrations is available as accessories.

### Specifications

- Basics of mechanical vibration, natural damped and forced vibrations
- Bar-type oscillator
- Three helical springs
- Unbalance exciter with DC motor, 0,77kg
- Displacement exciter with DC motor
- Electronic control unit with digital display, exciter frequency adjustable
- Oil-filled damper
- Electrically driven drum recorder
- Amplitude meter with electric contact for triggering equipment

### Technical Specifications

- Bar-type oscillator: LxWxH: 700x25x12mm, 1,6kg
- Helical springs
  - 0,75N/mm
  - 1,5N/mm
  - 3,0N/mm
- Exciter frequency: 0...50Hz, electronically controlled
- Unbalance of the unbalance exciter: 0...1000mmg
- Stroke of the displacement exciter: 20mm
- Damper constant: 5...15Ns/m, oil-filled
- Mechanical chart recorder

- Feed: 20mm/s
- Paper width: 100mm

### **Experiment Possibilities**

- Natural vibration
- Damped vibration
- Inertia force and displacement excitation
- Forced vibration
- Resonance
- Amplitude and phase response

### **Scope of Delivery**

- 1 frame
- 1 bar-type oscillator
- 3 helical springs
- 1 unbalance exciter
- 1 displacement exciter
- 1 exciter control unit
- 1 damper
- 1 amplitude meter
- 1 drum recorder
- 1 storage system with foam inlay
- 1 manual

Note: Specifications are subject to change.

### **Tesca Technologies Pvt. Ltd.**

IT-2013, Ramchandrapura Industrial Area, Sitapura Extension,  
Near Bombay Hospital, Vidhani Circle, Jaipur-302022, Rajasthan, India,  
Tel: +91-141-2771791 / 2771792; Email: info@tesca.in, tesca.technologies@gmail.com  
Website: www.tescaglobal.com

