


Hall Effect enables the charge carrier concentration and mobility to be determined by experiment. Direction of the Hall Voltage in silver indicates negative charge carriers, which is in agreement with concepts of the model of the 'free electron gas'. Limitations of this model are shown by the so called 'abnormal Hall Effect' of tungsten. The experiment carried out under identical conditions for tungsten show the Hall Voltage to have about same magnitude but opposite direction as in silver.

This can be explained by the 'Energy Band diagram'. The tungsten atom has $5s^2 5p^6 5d^4 6s^2$ electronic structure. When the atoms come close together to form the solid, the close lying states $5d$ and $6s$ broaden into bands, with s band broadening considerably more than the d band. This is because of the larger size of the s orbital. The figure schematically shows the allowed energies as a function of the interatomic distance. The number of allowed states is ten per atom in the d band and two in the s band. In tungsten



there are six electrons to be shared between these two bands. The result is that at the interatomic distance in tungsten there are holes in the d band and electrons in the s band, making tungsten predominantly a hole conductor.

This sort of mixed (electrons and holes) conduction is a general characteristic of transition metals. The apparatus consists of the following:

a) Hall Probes-Silver(HP-Ag)

Technical Specification

Material	: Silver(8X6X0.05 mm)
Contacts	: Press Type for current Spring type for Voltage
Hall Voltage	: ~17 mV/10A/10KG



b) Hall Probe-Tungsten (HP-W)

Technical Specification

Material	: Tungsten Strip (8 x 6 x 0.05 mm)
Contacts	: Press type for current Spring Type for Voltage
Hall Voltage	: ~15 mV/10A/10KG

Note: Specifications are subject to change.

Tesca Technologies Pvt. Ltd.

IT-2013, Ramchandrapura Industrial Area, Sitapura Extension,
Near Bombay Hospital, Vidhani Circle, Jaipur-302022, Rajasthan, India,
Tel: +91-141-2771791 / 2771792; Email: info@tesca.in, tesca.technologies@gmail.com
Website: www.tescaglobal.com

c) High Current Power Supply

Technical Specification

Range	: 0-20A continuously variable
Accuracy	: $\pm 0.5\%$
Regulation	: $\pm 0.5\%$ for $\pm 10\%$ variation of mains
Display	: 3½ digit, 7 Segment LED

d) Digital Microvoltmeter

It is a very versatile multipurpose instrument for the measurement of low dc voltage. It has 5 decade ranges from 1mV to 10mv with 100% over-ranging. For better accuracy and convenience, readings are directly obtained on 3½ digit DPM (Digital Panel Meter).

Technical Specifications

Range	: 1mV, 10mV, 100mV, 1V & 10V with 100% over-ranging.
Resolution	: 1µV
Accuracy	: $\pm 0.2\% \pm 1$ digit
Stability	: Within ± 1 digit
Input Impedance	: >1000MW (10MW on 10V range)
Display	: 3½ digit, 7 segment LED with autopolarity and decimal indication

e) Electromagnet



Technical Specification

Pole Pieces	: 75mm tapered to 25mm
Mag. Field	: 17KG $\pm 5\%$ at 10mm airgap
Energising Coils	: Two of approx. 13 each
Power	: 0-90Vdc, 3A, for coils in series 0-45Vdc, 6A, for coils in parallel

f) Constant Current Power Supply

g) Gaussmeter

The experiment is complete in all respect.



Note: Specifications are subject to change.

Tesca Technologies Pvt. Ltd.

IT-2013, Ramchandrapura Industrial Area, Sitapura Extension,
Near Bombay Hospital, Vidhani Circle, Jaipur-302022, Rajasthan, India,
Tel: +91-141-2771791 / 2771792; Email: info@tesca.in, tesca.technologies@gmail.com
Website: www.tescaglobal.com