

**Application:**

- Victor Meyer's apparatus is used to measure the vapour density of volatile organic compounds. A known mass of the compound is vaporised in an instrument called Victor Meyer tube. Vapours are obtained from the sample and it displaces an equal volume of air into a graduated tube.

Note: Specifications are subject to change, Photos shown above are Indicative, Actual Product can Vary.



Export Sales: +91-9829132777
India Sales: +91-9588842361



IT-2013, Ramchandrapura Industrial Area,
Sitapura Extension, Jaipur-302022, India.



info@tesca.in
www.tescaglobal.com