



55804 Experimental Set-Up has been designed specifically verification of Hook's Law. The set-up is complete in all respects and requires no other apparatus. Practical experience on this set up carries great educative value for Science and Engineering Students.

OBJECT

To verify the Hook's Law.

FEATURES

The Set-up consists of the followings :

- 01 Hook's Law Appratus :
Consisting of 12cm long mirror scale, spring, slotted weights with hanger & pointer, with complete stand.
- 02 Strongly supported by detailed Operating Instructions, giving details of Object, Theory, Design procedures, Report Suggestions and Book References.

OTHER APPARATUS REQUIRED

- 01 DIGITALSTOP CLOCK :
With START/STOP operation by means of toggle switch & RESET by a push button switch. It has a range of 999.9 seconds with resolution of 0.1 seconds and accuracy of $\pm 0.01\%$ (Quartz controlled). Display is thorough 4 no's of 12.5mm bright Seven Segment Displays and working voltage of the unit is $230V \pm 10\%$ 50Hz.

Note: Specifications are subject to change.

Tesca Technologies Pvt. Ltd.

IT-2013, Ramchandrapura Industrial Area, Sitapura Extension,
Near Bombay Hospital, Vidhani Circle, Jaipur-302022, Rajasthan, India,
Tel: +91-141-2771791 / 2771792; Email: info@tesca.in, tesca.technologies@gmail.com
Website: www.tesca.in

