



EEG 10-20 Method single Amplifier with USB Trainer Kit

- Real Time EEG Amplifier output, 10-20 Electrodes placements
- Electrodes placement select by Mode & Right & Left switch
- Study of Brain activity & signal processing
- EEG Left, Right & Self Test Mode, Gain Adjustable,
- Disc Electrodes Bipolar one electrodes set. Disc & Disposable Electrodes
- EEG DSO & USB output through 4mm socket & USB
- Software BLOWAVE, 6 Channel, Save, Open, Print, Excel Values, Measurement Voltage, Frequency, Time, Offset, Select 1 or 6 Channel, etc
- Test Point on Top of Instrument, 2mm Socket.
- Power Supply 230 \pm 10%, 50Hz, Dimensions 305 X100 X 250 mm \pm

Other Models

SN.	Order Code	Name of Item
01.	BM2121	EEG Amp. 3 Lead Electrodes
02.	BM2122	EEG Amp. 10-20 method Single 1Ch.
03.	BM2123	EEG Amp. 10-20 with Simulator 1Ch
04.	BM2124	EEG Amp. 10-20 with USB 1Ch
05.	BM2521	EEG Simulator

Operating Manual with working of instrument, operating & installation procedure, Block diagrams, connection diagrams, Electrodes placement, Standard signals & actual output signal, Calculations, Tables, Trouble shooting, Experiments list & etc.

Note: Specifications are subject to change.