



Output on DSO

## **EMG Biofeedback 3 Lead**

- Bipolar Real Time EMG Biofeedback 3 Lead Trainer
- No of Channel One, Also available 2, 3, 4 ... 8 Channels
- · Biofeedback available with USB connectivity
- EMG & Self Test Mode, Gain, Threshold & Audio adjustable
- Indication Mode Dot-Dot, Bar-Dot, Dot-Bar, Bar-Bar
- Disc or disposable Electrodes with Velcro.DSO output through 4mm socket
- Audio output CW 1wat ±. Test Point on Top of Instrument
- Power Supply 230 ±10%, 50Hz, Dimensions 305 X100 X 250 mm±

## Study on EMG Biofeedback

Study of Muscle activity, Muscle Response, with & Without Stimulator Responses & signal processing, Placement of Electrodes, Amplification, Filters, Gain, Self test, Voltage & Frequency responses & etc

## **Other Models**

| SN. | <b>Order Code</b> | Name of Item                    |
|-----|-------------------|---------------------------------|
| 01. | BM2131            | EMG Amp. 3 Lead with Audio      |
| 02. | BM2532            | EMG Simulator 3Lead             |
| 03. | BM2132            | EMG Amp. 3 Lead with Simulator  |
| 04. | BM2133            | EMG Amp. 3 Lead with Stimulator |
| 05. | BM2134            | EMG Biofeedback 3 Lead          |
| 06. | BM2135            | EMG Amp. with USB + Software    |
| 07. | BM2136            | EMG Amp. 2 channel              |
| 08. | BM2137            | EMG 3,4 Amplifiers (Channels)   |
| 09. | BM2138            | EGG Amp. 3 Lead with Audio      |

Operating Manual with working of instrument, operating & installation procedure, Block diagrams, connection diagrams, Electrodes placement, Standard signals & actual output signal, Calculations, Tables, Trouble shooting, Experiments list & etc.

Note: Specifications are subject to change.

Near Bombay Hospital, Vidhani Circle, Jaipur-302022, Rajasthan, India,

Tel: +91-9829132777; Email: info@tesca.in, tesca.technologies@gmail.com

© Website: www.tescaglobal.com

