



EEG Simulator Trainer Kit

Testing & Calibration of EEG Machine, Study of Signal Patterns.
Normal rhythms, Adjustable rhythms,
Rhythms Selection, Delta (d), Theta (f), Alpha (a), Beta (b), Gamma (g),
Amplitude adjustable 0 to 10mV±,
Output Lead A, G & R, Differential Output
Output socket 2mm or 4mm or Button Type,
Output to amplifier, recorder or EEG monitor or DSO
Red, Blue & White light for Eyes activity
Power Supply 250V, 50Hz, Dimensions 250 X 75 X 250mm±

Other EMG Models

SN.	Order Code	Name of Item
01.	BM2121	EEG Amp. 3 Lead Electrodes
02.	BM2122	EEG Amp. 10-20 method Single 1Ch.
03.	BM2123	EEG Amp. 10-20 with Simulator 1Ch
04.	BM2124	EEG Amp. 10-20 with USB 1Ch
05.	BM2521	EEG Simulator

Operating Manual with working of instrument, operating & installation procedure, Block diagrams, connection diagrams, Electrodes placement, Standard signals & actual output signal, Calculations, Tables, Trouble shooting, Experiments list,

Note: Specifications are subject to change.