



Applications:

Used to determine the particle charge of Bitumen Emulsion.

Apparatus Details:

- The equipment comprises of following items
- Digital Milliammeter scale up to 10 mA on support base
- · Variable resistor
- Two stainless steel electrodes
- · Insulating device
- Glass rod
- Rehostat-2000 ohms
- Dimensions : 200 x 200 x 600mm

Technical specification

• Overall Dimensions: 700 X 350 X 240(L X W X D) mm

: 250 ml Beaker Battery : 12 v

Ammeter : 0.1 amp
Power Supply : 250 Volts, 1ph, 50/60Hz
Weight : 6 kg approximately

Note: Specifications are subject to change.

Near Bombay Hospital, Vidhani Circle, Jaipur-302022, Rajasthan, India, Tel: +91-9829132777; Email: info@tesca.in, tesca.technologies@gmail. Tel: +91-9829132777; Email: info@tesca.in, tesca.technologies@gmail.com