



## Description :

- 50 Strips To Screen Nitrate, pH, SG, Blood, Glu, Bil, Urobilinogen, Ketone, leucocytes, Protein in Urine Packed In Airtight Bottle And Box Protection With Silica. Good Till Opened Kindly read the instructions properly before use, you should not keep the strip dipped in urine for more than 1/2 minute and afterwards you have to lay the strip on the table horizontally and ensure there is no much sample coming out of strips, you do have to takeout excess sample from the strip without touching it by tapping the strip on the container or on washbasin

## Technical Specifications :

- pH Level : 5-9
- Net Quantity : 50 strips

*Note: Specifications are subject to change, Photos shown above are Indicative, Actual Product can Vary.*



Export Sales: +91-9829132777  
India Sales: +91-9588842361



IT-2013, Ramchandrapura Industrial Area,  
Sitapura Extension, Jaipur-302022, India.



info@tesca.in  
www.tescaglobal.com