

**Specifications:**

Complete demonstration of Pascal's Law. The apparatus shows that a change in pressure in fluids is transmitted equally in all directions to all parts of a fluid. The device contains a hollow sphere connected to a hand pump with plunger that applies pressure to fluid in the ball. The fluid exits all 10 equidistant nozzles in the ball with equal force in all directions. With instructions.

Note: Specifications are subject to change, Photos shown above are Indicative, Actual Product can Vary.



Export Sales: +91-9829132777
India Sales: +91-9588842361



IT-2013, Ramchandrapura Industrial Area,
Sitapura Extension, Jaipur-302022, India.



info@tesca.in
www.tescaglobal.com