



Order Code : 20213547.5.6

Name : Milkan oil drop apparatus



Specifications:

- Millikan's oil drop experiment measured the charge of an electron. Before this experiment, existence of subatomic particles was not universally accepted. Millikan's apparatus contained an electric field created between a parallel pair of metal plates, which were held apart by insulating material.

Note: Specifications are subject to change, Photos shown above are Indicative, Actual Product can Vary.

TESCA TECHNOLOGIES PVT. LTD.

IT-2013, Ramchandrapura Industrial Area, Sitapura Extension, Jaipur-302029, Rajasthan, India.
Ph/ Fax: 91-141-2771791, 2771792; Email: info@tesca.in, tesca.technologies@gmail.com
Website: www.tesca.in