



DESCRIPTION:

The apparatus consists of a glass tube with one end having bell mouth entrance; connected to a constant head water tank, at the other end a valve is provided to vary the flow rate. The tank is of sufficient capacity to store water. A capillary tube is introduced centrally in the bell mouth for feeding dye from a small container placed at the top of tank, through polythene tubing. By varying the rate of flow, the Reynold's number is changed. This also changes the type of flow. Visual observation of dye (Thread) will indicate the type of flow, which can be confirmed from the Reynold's number computed. Present set-up is self-contained water re-circulating unit, provided with a sump tank and a centrifugal pump etc. Flow control valve and by-pass valve are fitted in water line. Flow rate of water is measured with the help of measuring cylinder and stop watch.

Sensor for measuring discharge of water. Data acquisition system with USB for connection to PC

Technical Details:

- Tube: Material Borosilicate Glass
- Dye vessel: Material Stainless Steel, Suitable Capacity
- Capillary Tube: Material Copper/Stainless Steel.
- Constant Head Water Tank: Capacity 40 Ltrs.
- Water Circulation: FHP Pump.
- Flow Measurement: Using Measuring Cylinder 2 LITERS
- Flow sensor for measuring discharge
- Sump Tank: Capacity 60 Ltrs.
- Stop Watch: Electronic.
- Control Panel Comprises of : Standard make On/Off Switch, Mains Indicator, etc.
- Instruction Manual: An ENGLISH instruction manual will be provided along with the Apparatus.
- Tanks will be made of Stainless Steel.
- The whole set-up is well designed and arranged in a good quality painted structure.

Experiment:

- To determine the Reynold's number and hence the type of flow either laminar or turbulent.

Note: Specifications are subject to change, Photos shown above are Indicative, Actual Product can Vary.



Export Sales: +91-9829132777
 India Sales: +91-9588842361



IT-2013, Ramchandrapura Industrial Area,
 Sitapura Extension, Jaipur-302022, India.



info@tesca.in
 www.tescaglobal.com

- To study transition zone.

Utilities Required:

- Electricity Supply: Single Phase 220 V AC, 50Hz, 5-15 amp combined socket with earth connection.
- Water Supply (Initial Fill).
- Floor Drain required.
- Floor Area Required: 1.5 m x 0.75 m
- Chemical Required: Dye (KMnO₄) - 10 gm

Note: Specifications are subject to change, Photos shown above are Indicative, Actual Product can Vary.



Export Sales: +91-9829132777
India Sales: +91-9588842361



IT-2013, Ramchandrapura Industrial Area,
Sitapura Extension, Jaipur-302022, India.



info@tesca.in
www.tescaglobal.com