



Order Code : 20213547.4.5

Name : Modern optics apparatus



Specifications:

Fresnel Biprism setup illustrates about the phenomenon of interference of light. With the help of interference phenomenon, we can find the wavelength of monochromatic light source like sodium light in a very easy manner. Complete with optical bench of one meter long

Note: Specifications are subject to change, Photos shown above are Indicative, Actual Product can Vary.

TESCA TECHNOLOGIES PVT. LTD.

IT-2013, Ramchandrapura Industrial Area, Sitapura Extension, Jaipur-302029, Rajasthan, India.

Ph/ Fax: 91-141-2771791, 2771792; Email: info@tesca.in, tesca.technologies@gmail.com

Website: www.tesca.in