



Technical Specifications :

- Nigrosin Stain, 10% w/v, is supplied by Micromaster. It is available in the form of a blackish violet coloured solution which is clear without any particles. It is advisable to store this product between 10- 30C in tightly closed container and away from bright light.
- Quantity provided: 100ml/bottle
- Nigrosin stain, 10% w/v is used as staining solution for negative staining. Negative staining method permits visualization of the usually transparent and unstainable capsule of many organisms, most importantly *Cryptococcus neoformans*. Nigrosin consists of a suspension of fine particles of carbon. These form a dark background, against which capsules are clearly seen as a result of displacement of the carbon particles.

Note: Specifications are subject to change, Photos shown above are Indicative, Actual Product can Vary.



Export Sales: +91-9829132777
India Sales: +91-9588842361



IT-2013, Ramchandrapura Industrial Area,
Sitapura Extension, Jaipur-302022, India.



info@tesca.in
www.tescaglobal.com