



EOG Amplifier 3 Lead with USB BIOWAVE

- Bipolar Real Time EOG Amplifier 3 Lead Trainer
- 3 Lead Electrodes placements Active, Reference & Ground. Electrodes Disc with Velcro & Disposable Electrodes. Mode Self Test & EOG, Gain & Self-Test Adjustable. LED indication, DSO output through 4mm socket. Test Point on Top.
- Power Supply 230 \pm 10%, 50Hz, Dimensions 305 X100 X 250 mm \pm

Study on EOG Amplifier

Study of electric potentials generated by the eyes movements of a stationary subject. Detection Moment Right, Left, Up & Down. Study of Pick voltage & Duration of Moments & signal processing Placement of Electrodes, Amplification, Filters, Gain, Self test, Voltage & Frequency responses & etc

BIOWAVE Software

Software for measurement of EOG single. To measure Voltage, Time & Frequency of Signal Pattern. Setting Voltage & Time Expand the signal, Save As, Open, Print, etc. For advance study Excel values, offset setting, etc

Other Models

SN.	Order Code	Name of Item
01.	BM2126	EOG Amplifier 3 Lead Electrodes
02.	BM2127	EOG Amplifier 3 Lead with Simulator
03.	BM2128	EOG Amplifier 3 Lead with USB BIOWAVE
04.	BM2526S	EOG Simulators
05.	BM2527	EOG Simulator with light simulator

Operating Manual with working of instrument, operating & installation procedure, Block diagrams, connection diagrams, Electrodes placement, Standard signals & actual output signal, Calculations, Tables, Trouble shooting, Experiments list & etc.

Note: Specifications are subject to change.