



Applications:

To prepare test specimen by molding in a specially designed sample die.

Details:

- Manually operated hydraulic press of 5000 kg (50 KN) capacity
- Pressure applied on the sample with help of manually operated Hydraulic Jack
- Platen sizes 15 x 15 cms approximately
- Supplied complete with sample die of 50mm (Optionally other dies can be supplied at extra cost)

Note: ELECTRICALLY OPERATED BUTTON PRESSURE UP TO 3000 KN CAPACITY ALSO AVAILABLE.wheels, will be supplied for bigger version autoclave. It will cost extra. The price on request.

Note: Specifications are subject to change.

Canal Strand Str

Near Bombay Hospital, Vidhani Circle, Jaipur-302022, Rajasthan, India,

Tel: +91-9829132777; Email: info@tesca.in, tesca.technologies@gmail.com

[∞] Tel: +91-9029192,..., _ [∞] Website: www.tescaglobal.com