



Features:-

Idea of ready to use, convenient, cost effective & easy to handle test strips to detect the chemical adulterants & establish the Microbial Quality (freshness) of milk. Handy tool kit to detect harmful adulterants and assess quality of milk consumed daily. Test can be conducted in 3 easy steps – Pour milk, insert test strip and observe change in color. Testing does not require any electricity and add on support.

Note: Specifications are subject to change, Photos shown above are Indicative, Actual Product can Vary.



Export Sales: +91-9829132777
India Sales: +91-9588842361



IT-2013, Ramchandrapura Industrial Area,
Sitapura Extension, Jaipur-302022, India.



info@tesca.in
www.tescaglobal.com